

## **Double staining of lymphatic endothelial cells and HUVECs**

*(Protocol for cell culture dish)*

### **Material**

- PBS<sup>++</sup>
- 4% PFA in PBS
- 0.1% Triton/BSA in PBS
- 1% BSA in PBS<sup>++</sup>
- DAPI/PBS (1µg/ml)

### **Antibodies** (for example!)

- anti-human Prox-1, Rabbit, use 1:1000 (RELIATech Cat# 102-PA32)
- anti-human CD31

### **Protocol**

- wash cells 2x with PBS
- fix cells for 5min at RT with 4% PFA in PBS
- wash cells 2x with PBS
- incubate for 45min in 1% BSA/PBS at RT
- perforation of membrane by incubation in 0.1% Triton/BSA for 5min
- wash cells 2x with PBS
- incubation with primary antibody for 2h at RT
- wash cells 2x with PBS
- incubation with secondary antibody for 1h at RT
- wash cells 2x with PBS
- incubation for 10min at RT in DAPI (optional)
- cover cells with Flouromount directly in the culture dish

☞ samples are ready for analysis with the microscope!